

**TOP
LONGMAX**[®]

**A RELIABLE SOURCE
OF WHEELCHAIRS &
RELATED PRODUCTS**



COMPANY PROFILE

TOP LONGMAX was established in 1984. Specialize in producing medical and orthopedic appliances, such as Wheelchair ; Folding Chair ; Walker ; Shower Chair ; Rollator ; Commode ; Quad Cane ; Crutch and all sorts of Support and Other related optional parts.

Our company creates new products constantly, hope to provide the most comfortable and fit in with humanity equipment. Strict quality control allows us to provide products to stringent markets worldwide. We are proud of our excellent products quality and superior services which have approved from our customers.

We had plentiful technique and experience to match with you requires.

PRODUCTS INFORMATION

Wheelchair	Shower Chair
C. P. Chair	Grab Bar
Power Chair	Overbed Table
Rollator	Stretcher
Commode	Trolley
Walker	Safety Belt
Quad Cane	Medical equipments
Crutch	
Cushion	Support ... New Catalogue
Bath Bench	Walking Cane ... New Catalogue

WHEELCHAIR CODE INFORMATION

TOP - **CASTER TYPE** **REAR WHEEL TYPE** - **ARM TYPE** **FOOTREST TYPE** - **CHAIR SIZE**

CASTER TYPE

S : 8" x 1" **Solid** caster

P : 8" x 2" **Pneumatic** caster

U : 8" x 2" **PU** flatless caster

REAR WHEEL TYPE

S : 24" x 1" **Solid** rear wheel

P : 24" x 1-3/8" **Pneumatic** rear wheel

U : 24" x 1-3/8" **PU** flatless rear wheel

ARM TYPE

10 : Detachable **Desk** armrest

15 : Detachable **Waterfall Desk** armrest

20 : **Fixed** full length armrest

30 : Detachable **Full Length** armrest

35 : Detachable **Waterfall Full Length** armrest

FOOTREST TYPE

AF : Adjustable height **Fixed** footrest

SDF : **Swingaway Detachable Footrest**

SDEL : **Swingaway Detachable Elevating Legrest**

CHAIR SIZE

A : Adult size **18"**

N : Narrow size **16"**

C : Child size **14"**

HD : Heavy Duty size **20"**

TOP-PP-35SDEL-AD

Deluxe Style Wheelchair

- ▲ Detachable waterfall full length armrest
- ▲ Swing away and detachable elevating legrest
- ▲ Deluxe upholstery seat and back
- ▲ 8" x 2" pneumatic front caster
- ▲ 24" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDF / 10 ; 15



TOP-SS-15SDEL-AD

Deluxe Style Wheelchair

- ▲ Detachable waterfall desk armrest
- ▲ Swing away and detachable elevating legrest
- ▲ Deluxe upholstery seat and back
- ▲ 8" x 1" solid front caster
- ▲ 24" x 1-3/8" solid rear wheel
- ▲ Heavy duty brake

Options

- PP ; SP ; PS ; US ; UP ; UU / SDF / 10 ; 30 ; 35



TOP-PP-35SDF-AD

Deluxe Style Wheelchair

- ▲ Detachable waterfall full length armrest
- ▲ Swing away and detachable footrest
- ▲ Deluxe upholstery seat and back
- ▲ 8" x 2" pneumatic front caster
- ▲ 24" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDEL / 10 ; 15 ; 35



TOP-PP-15SDF-AD

Deluxe Style Wheelchair

- ▲ Detachable waterfall desk armrest
- ▲ Swing away and detachable footrest
- ▲ Deluxe upholstery seat and back
- ▲ 8" x 2" pneumatic front caster
- ▲ 24" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDEL / 10 ; 30 ; 35



TOP-PP-30SDEL-A

Universal Style Wheelchair

- ▲ Detachable full length armrest
- ▲ Swing away and detachable elevating legrest
- ▲ PVC upholstery seat and back
- ▲ 8" x 2" pneumatic front caster
- ▲ 24" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDF / 10



TOP-SS-10SDEL-A

Universal Style Wheelchair

- ▲ Detachable desk armrest
- ▲ Swing away and detachable elevating legrest
- ▲ PVC upholstery seat and back
- ▲ 8" x 1" solid front caster
- ▲ 24" x 1-3/8" solid rear wheel
- ▲ Heavy duty brake

Options

- PP ; SP ; PS ; US ; UP ; UU / SDF / 30



TOP-PP-30SDF-A

Universal Style Wheelchair

- ▲ Detachable full length armrest
- ▲ Swing away and detachable footrest
- ▲ PVC upholstery seat and back
- ▲ 8" x 2" pneumatic front caster
- ▲ 24" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDEL / 10



TOP-PP-10SDF-A

Universal Style Wheelchair

- ▲ Detachable desk armrest
- ▲ Swing away and detachable footrest
- ▲ PVC upholstery seat and back
- ▲ 8" x 2" pneumatic front caster
- ▲ 24" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDEL / 30



TOP-SS-20SDF-A **Standard Style Wheelchair**

- ▲ Fixed full length armrest
- ▲ Swing away and detachable footrest
- ▲ PVC upholstery seat and back
- ▲ 8" x 1" solid front caster
- ▲ 24" x 1-3/8" solid rear wheel
- ▲ Standard brake

Options

- PP ; SP ; PS ; US ; UP ; UU / SDEL



TOP-SS-20AF-A **Economy Style Wheelchair**

- ▲ Fixed full length armrest
- ▲ Adjustable height fixed footrest
- ▲ PVC upholstery seat and back
- ▲ 8" x 1" solid front caster
- ▲ 24" x 1-3/8" solid rear wheel
- ▲ Standard brake

Options

- PP ; SP ; PS ; US ; UP ; UU



TOP-PP-01A **Folding Transport Chair**

- ▲ Detachable desk armrest
- ▲ Swing away and detachable footrest
- ▲ PVC upholstery seat and back
- ▲ 8" x 2" pneumatic front caster
- ▲ 12" x 2-1/2" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDEL / 30



TOP-PP-30SDF-AP **Light Weight Wheelchair**

- ▲ Flip-up armrest
- ▲ Swing away and detachable footrest
- ▲ Nylon upholstery seat and back
- ▲ 7" x 1-3/4" pneumatic front caster
- ▲ 24" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDEL



TOP-PP-10SDEL-AR

Full Reclining Wheelchair

- ▲ Detachable desk armrest
- ▲ Swing away and detachable elevating legrest
- ▲ Full 90° reclining backrest
- ▲ Removable headrest
- ▲ 8" x 2" pneumatic front caster
- ▲ 24" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU / SDF / 30



TOP-PP-20AF-AMV

Folding Arm Moving Chair

- ▲ Fixed armrest
- ▲ Adjustable height fixed footrest
- ▲ Two arms moving, forward, backward, turn right and turn left smoothly
- ▲ 8" x 2" pneumatic front caster
- ▲ 20" x 1-3/8" pneumatic rear wheel
- ▲ Heavy duty brake



TOP-SS-20AF-AC

Economy Folding Commode Wheelchair

- ▲ Fixed armrest
- ▲ Adjustable height fixed footrest
- ▲ Removable headrest
- ▲ Partial board toilet seat with cushion and removable stainless steel pail
- ▲ 8" x 1" solid front caster
- ▲ 24" x 1-3/8" solid rear wheel
- ▲ Heavy duty brake

Options

- SS ; SP ; PS ; US ; UP ; UU



Attendant table



T-301

Made by Wooden



T-302

Made by Acrylic

TOP-AL-30SDF-A1

Aluminum Alloy Wheelchair

- ▲ Detachable and adjustable foam padded armrest
- ▲ Swingaway and detachable footrest
- ▲ Nylon upholstery seat and back, with seat foam cushion
- ▲ Detachable and adjustable height backrest
- ▲ 8" x 1" pneumatic front caster
- ▲ 24" x 1-3/8" quick release pneumatic rear wheel
- ▲ Adjustable front and rear wheel base
- ▲ Adjustable heavy duty brake



TOP-AL-30SDF-A1A

Aluminum Alloy Wheelchair

- ▲ Detachable and adjustable foam padded armrest
- ▲ Swing away and detachable footrest
- ▲ Nylon upholstery seat and back, with seat foam cushion
- ▲ Detachable and adjustable height backrest
- ▲ 8" x 1" PU front caster
- ▲ 24" x 1-3/8" quick release pneumatic rear wheel with plastic handrim
- ▲ Adjustable heavy duty brake



TOP-AL-30SDF-A2

Aluminum Alloy Wheelchair

- ▲ Detachable flip backward armrest
- ▲ Swingaway and detachable footrest
- ▲ Nylon upholstery seat and back, with seat foam cushion
- ▲ Detachable and adjustable height backrest
- ▲ 8" x 1" pneumatic front caster
- ▲ 24" x 1-3/8" quick release pneumatic rear wheel
- ▲ Adjustable front and rear wheel base
- ▲ Adjustable heavy duty brake



TOP-AL-30SDF-A2A

Aluminum Alloy Wheelchair

- ▲ Detachable flip backward armrest
- ▲ Swing away and detachable footrest
- ▲ Nylon upholstery seat and back, with seat foam cushion
- ▲ Detachable and adjustable height backrest
- ▲ 8" x 1" PU front caster
- ▲ 24" x 1-3/8" quick release pneumatic rear wheel
- ▲ Mag. rear wheel with steel handrim
- ▲ Adjustable heavy duty brake



TOP-AL-30SDF-A5

Aluminum Alloy Wheelchair

- ▲ Detachable flip backward armrest
- ▲ Swingaway and detachable footrest
- ▲ Nylon upholstery seat and back, with seat foam cushion
- ▲ Detachable and adjustable height backrest
- ▲ 8" x 1" pneumatic front caster
- ▲ 24" x 1-3/8" quick release pneumatic rear wheel
- ▲ Adjustable heavy duty brake



TOP-AL-20AFJ-AF

Aluminum Alloy Wheelchair

- ▲ Fixed armrest
 - ▲ Adjustable height fixed footrest
 - ▲ Nylon upholstery seat and back
 - ▲ Foldable backrest
 - ▲ 6" PU front caster
 - ▲ 22" spoked pneumatic rear wheel
- Options
- Attendant control brake lever (TOP-AL-20AFJ-AFH)
 - Fix backrest



TOP-AL-01FB-14"

Aluminum Transporter Wheelchair

- ▲ Fixed armrest
 - ▲ Adjustabel height fixed footrest
 - ▲ Deluxe upholstery seat and back
 - ▲ Seat width 14"
 - ▲ 6" solid front caster
 - ▲ 16" pneumatic rear wheel
- Options
- Fixed backrest (TOP-AL-01AB-14")



TOP-AL-20AF-16"

Economy Aluminum Alloy Wheelchair

- ▲ Fixed armrest
 - ▲ Adjustabel height fixed footrest
 - ▲ Deluxe upholstery seat and back
 - ▲ Seat width 16"
 - ▲ 8" x 2" PU front caster
 - ▲ 24" x 1-3/8" pneumatic rear wheel with plastic handrim
- Options
- Foldable backrest (TOP-AL-20AF-16"A)
 - Seat width 18"



TOP-AL-SPORT

Aluminum Alloy Sport Wheelchair

- ▲ Fixed frame
- ▲ Quick knock down and assembly
- ▲ Adjustable and fold down backrest
- ▲ 4" quick release PU front caster
- ▲ 24" x 1" quick release high pressure rear wheel
- ▲ Heavy duty aluminum alloy brake



TOP-AL-SPORT-A

Aluminum Alloy Sport Wheelchair

- ▲ Super light 8.90kgs
- ▲ Adjustable seat height 44~52cm from floor
- ▲ 4" quick release PU front caster
- ▲ 24" x 1" quick release high pressure rear wheel
- ▲ Heavy duty aluminum alloy brake



TOP-02D

Deluxe Type Commode

- ▲ Swing down armrest
- ▲ Swing away and detachable footrest
- ▲ Molded toilet seat with cushion
- ▲ 5" PU front caster
- ▲ 5" PU rear wheel with brake



TOP-AL-04

Aluminum Alloy Commode / Shower Chair

- ▲ Adjustable height and movable armrest
 - ▲ Swing away and detachable footrest
 - ▲ PVC upholstery back & PU molded seat
 - ▲ 5" PU front caster
 - ▲ 22" x 1-3/8" PU rear wheel
 - ▲ Heavy duty aluminum alloy brake
- Options
- 5" PU rear wheel with brakes (TOP-AL-04B)



TOP-AL-04B





TOP-747 Alloy Aluminum Aircraft Wheelchair

- ▲ Lightweight aircraft-grade aluminum alloy frame
- ▲ Sturdy non-foldable frame
- ▲ Flip-up armrest
- ▲ Adjustable height headrest
- ▲ Detachable and adjustable height backrest
- ▲ Four way safety belt, shank belt and thigh belt
- ▲ Seat width 13"
- ▲ 6" x 1" solid front caster
- ▲ 8" x 1" solid rear wheel
- ▲ Two 3" wheels under the footrest for more stably
- ▲ Heavy duty brakes on the rear



TOP-PP-10SDF-CH Folding Child Wheelchair

- ▲ Detachable desk armrest
- ▲ Swing away detachable and vertical adjustable height footrest
- ▲ 8" x 2" pneumatic front caster
- ▲ 22" Mag. rear wheel with plastic handrim
- ▲ Heavy duty brake

Options

- Seat width : 10' ; 12' ; 13' ; 14' ; 16'
- Chrome plated steel ; Powder coating
- 22"spoke rear wheel with metal handrim



TOP-505 Correction and Positioning Wheelchair

- ▲ Aluminum alloy frame
- ▲ Vertical adjustable height footrest
- ▲ Full 90° reclining backrest
- ▲ Adjustable anti-tipper and headrest
- ▲ 5" x 1" PU front caster
- ▲ 12-1/2" x 2-1/4" quick release PU rear wheel
- ▲ Safety special seat belt and leg strap

Options

- Seat width : 12" ; 13" ; 14" ; 15" ; 16"
- Folding backrest
- Attendant table

Attendant table





TOP-506

Correction and Positioning Wheelchair

- ▲ Aluminum alloy frame
- ▲ Swing away detachable and adjustable height footrest
- ▲ Full 90° reclining and folding backrest
- ▲ Adjustable anti-tipper and headrest
- ▲ 5" x 1" PU front caster
- ▲ 12-1/2" x 2-1/4" quick release PU rear wheel
- ▲ Safety special seat belt and heel strap

Options

- Seat width : 12" ; 13" ; 14" ; 15" ; 16"
- Fixed backrest
- 22" x 1-3/8" PU rear wheels (TOP-506A)
- Attendant table



TOP-506B

Correction and Positioning Wheelchair

- ▲ Aluminum alloy frame
- ▲ Swing away detachable and adjustable height footrest
- ▲ Full 90° reclining and folding backrest
- ▲ Adjustable armrest and headrest
- ▲ Hard seat base and backrest with memory soft cushion
- ▲ 5" x 1" PU front caster
- ▲ 12-1/2" x 2-1/4" quick release PU rear wheel
- ▲ Safety special seat belt and leg strap

Options

- Seat width : 12" ; 13" ; 14" ; 15" ; 16"
- Fixed backrest
- 22" x 1-3/8" PU rear wheels (TOP-506BA)
- Attendant table



TOP-507A

Correction and Positioning Wheelchair

- ▲ Aluminum alloy frame
- ▲ Adjustable height armrest
- ▲ Vertical adjustable height footrest
- ▲ Full 90° reclining and folding backrest
- ▲ Adjustable anti-tipper and headrest
- ▲ 6" x 1-1/4" PU front caster
- ▲ 22" x 1-3/8" quick release PU rear wheel
- ▲ Safety special seat belt and leg strap
- ▲ Sparately function for seat and backrest reclined (Can be backrest reclined only or seat and backrest reclined together)

Options

- Seat width : 12" ; 13" ; 14" ; 15" ; 16"
- Fixed backrest
- 12-1/2" x 2-1/4" PU rear wheels (TOP-507)
- Attendant table

Attendant table





TOP-601

Reclining Stander

According to the motion developmental delay, provide full level to 90° standing training and fix the joint and knee at straighten position, improve patient head, body, knee, ankle...etc, get full protect also can avoid muscle and joint contracture.

- ▲ Dimension : (A) Height 75cm-120cm
(B) Height 120cm-160cm
- ▲ Bed : H 74cm x L 135 B160cm x W 35 cm
- ▲ Base : 96 x 72cm full level to 90°
- ▲ Footrest to railing height 82cm-105cm
- ▲ Footrest can be adjusted separately
- Options
 - Body band, pelvis band, knee band and leg band
 - Railing table



TOP-602

Vertical Stander

According to the cerebral palsy, dyskinesia patient to fix the joint and knee at straighten position, improve patient hand and body control force, avoid muscle and joint contracture.

- ▲ Dimension : (A) Height 75cm-120cm
(B) Height 120cm-160cm
- ▲ Adjustable angle : 70°-90°
- ▲ Removable table and body side support can change height and wide.
- ▲ Separately knee supporter full change side, depend on children age.
- Options
 - Body band, waist band, pelvis band, pelvis holder, feet band

TOP-603 Child Cart



- ▲ Seat : W 8"-14" ; D 9"-15"
- ▲ Chair : H 25" ; 90-135 degree
- ▲ Fold Size : 12" x 13" x 48"
- ▲ Footrest : 7"
- ▲ Front wheels are 6"
- ▲ Back wheels are 12"
- ▲ Loading capacity : 132lbs, 60kg

TOP-604 Stainless Steel Vertical Standing Support



Adjustable for height, the pedal can be adjusted separately for height, the pelvis band and knee joint fixing board can be fixed at appropriate position.

- ▲ Complete material is stainless steel
- ▲ Dimension : L 70cm x W 70cm x H 100cm-140cm
- ▲ Height of the support adjustable 17cm-57cm

TOP-POWER-A1

Folding Power Wheelchair



- ▲ Detachable and adjustable height armrest
- ▲ Swing away and detachable footrest
- ▲ Rigid one piece side frame made by powder coating steel
- ▲ Unique mechanism to easily mount and dismount batteries
- ▲ Programmable controller with wide range of selection
- ▲ Tire caster : 8" x 2" foam filled tire
- ▲ Drive rear wheel : 12-1/2" x 2-1/4" foam filled tire
- ▲ Capacity : 250 lbs
- ▲ Max speed : 6.40 km (4 MPH)
- ▲ Turning radius : 58 cm (23")
- ▲ Range : 32 km (20 miles)
- ▲ Controller : P&G V.S.L. 50505 or Dynamic DA50-C51
- ▲ Motor : DC 24V permanent magnetic motor, 250W 3800rpm
- ▲ Safe climbing angle : 12 degrees
- ▲ Brake : Intelligent, regenerative, electromagnet brake
- ▲ Batteries : 12V U1 (36AH) x 2
- ▲ Battery charger : 5 Amp off-board

Battery charger



TOP-POWER-A2

Folding Power Wheelchair



- ▲ Detachable and adjustable height armrest
- ▲ Swing away and detachable footrest
- ▲ Rigid one piece side frame made by powder coating steel
- ▲ Unique mechanism to easily mount and dismount batteries
- ▲ Programmable controller with wide range of selection
- ▲ Tire caster : 8" x 2" foam filled tire
- ▲ Drive rear wheel : 24" x 1-3/8" PU tire
- ▲ Capacity : 250 lbs
- ▲ Max speed : 10.50 km (7 MPH)
- ▲ Turning radius : 58 cm (23")
- ▲ Range : 32 km (20 miles)
- ▲ Controller : P&G V.S.L. 50505 or Dynamic DA50-C51
- ▲ Motor : DC 24V permanent magnetic motor, 250W 3800rpm
- ▲ Safe climbing angle : 12 degrees
- ▲ Brake : Intelligent, regenerative, electromagnet brake
- ▲ Batteries : 12V U1 (36AH) x 2
- ▲ Battery charger : 5 Amp off-board



TOP-POWER-A3

Suspension Power Wheelchair



- ▲ Detachable and adjustable height armrest
- ▲ Folding and adjustable angle backrest
- ▲ Rigid frame made by alloy aluminum
- ▲ Unique mechanism to easily mount and dismount batteries
- ▲ Programmable controller with wide range of selection
- ▲ Front and rear suspension for rough terrain
- ▲ Tire caster : 9" x 2.80"(2.5/2.8-4) foam filled tire
- ▲ Drive rear wheel : 16" x 4" (4.0-8) foam filled tire
- ▲ Capacity : 250 lbs
- ▲ Max speed : 8km (4 MPH)
- ▲ Turning radius : 58 cm (23")
- ▲ Range : 32 km (20 miles)
- ▲ Controller : P&G V.S.L. 50505 or Dynamic DA50-C51
- ▲ Motor : DC 24V permanent magnetic motor, 250W 3800rpm
- ▲ Safe climbing angle : 12 degrees
- ▲ Brake : Intelligent, regenerative, electromagnet brake
- ▲ Batteries : 12V U1 (36AH) x 2
- ▲ Battery charger : 5 Amp off-board

Options

- Safe bumper with suspension



TOP-POWER-A3H

Suspension Power Wheelchair



- ▲ Detachable and adjustable height armrest
- ▲ Folding and adjustable angle backrest
- ▲ Rigid frame made by alloy aluminum
- ▲ Unique mechanism to easily mount and dismount batteries
- ▲ Programmable controller with wide range of selection
- ▲ Front and rear suspension for rough terrain
- ▲ Tire caster : 9" x 2.80" (2.5/2.8-4) foam filled tire
- ▲ Drive rear wheel : 16" x 4" (4.0-8) foam filled tire
- ▲ Capacity : 250 lbs
- ▲ Max speed : 7.60 km (5 MPH)
- ▲ Turning radius : 58 cm (23")
- ▲ Range : 32 km (20 miles)
- ▲ Controller : P&G V.S.L. 50505 or Dynamic DA50-C51
- ▲ Motor : DC 24V permanent magnetic motor, 250W 3800rpm
- ▲ Safe climbing angel : 12 degrees
- ▲ Brake : Intelligent, regenerative, electromagnet brake
- ▲ Batteries : 12V U1 (36AH) x 2
- ▲ Battery charger : 5 Amp off-board

Options

- Safe bumper with suspension



TOP-POWER-AIR

Reclining Power Wheelchair



- ▲ Detachable and flip-up adjustable height armrest
- ▲ Comfort gel with memory foam seat and back cushion
- ▲ Folding and adjustable angle backrest
- ▲ Full 90° reclining backrest
- ▲ Detachable and adjustable headrest
- ▲ Rigid one piece side frame made by powder coating steel
- ▲ Unique mechanism to easily mount and dismount batteries
- ▲ Programmable controller with wide range of selection
- ▲ Tire caster : 8" x 2" foam filled tire
- ▲ Drive rear wheel : 12-1/2" x 2-1/4" foam filled tire
- ▲ Max speed : 6.50 km (4 MPH)
- ▲ Turning radius : 58 cm (23")
- ▲ Range : 32 km (20 miles)
- ▲ Controller : P&G V.S.L. 50505 or Dynamic DA50-C51
- ▲ Motor : DC 24V permanent magnetic motor, 320W 3500rpm
- ▲ Safe climbing angle : 12 degrees
- ▲ Brake : Intelligent, regenerative, electromagnet brake
- ▲ Batteries : 12V U1 (36AH) x 2
- ▲ Battery charger : 5 Amp off-board

Battery charger



TOP-POWER-A1C

Chin Control Power Wheelchair



- ▲ Detachable and flip-up adjustable height armrest
- ▲ Comfort gel with memory foam seat and back cushion
- ▲ Folding and adjustable angle backrest
- ▲ Adjustable chin control
- ▲ Detachable and adjustable headrest
- ▲ Rigid one piece side frame made by powder coating steel
- ▲ Unique mechanism to easily mount and dismount batteries
- ▲ Programmable controller with wide range of selection
- ▲ Tire caster : 8" x 2" foam filled tire
- ▲ Drive rear wheel : 12-1/2" x 2-1/4" foam filled tire
- ▲ Max speed : 6.40 km (4 MPH)
- ▲ Turning radius : 58 cm (23")
- ▲ Range : 32 km (20 miles)
- ▲ Controller : P&G or Dynamic
- ▲ Motor : DC 24V permanent magnetic motor, 250W 3800rpm
- ▲ Safe climbing angle : 12 degrees
- ▲ Brake : Intelligent, regenerative, electromagnet brake
- ▲ Batteries : 12V U1 (36AH) x 2
- ▲ Battery charger : 5 Amp off-board



Options

- Safe bumper with suspension



TOP-POWER-A4

Solar Power and Electric Two Way Use Wheelch



- ▲ Max power rating p_{max} (WP) : $55W \times 2 = 110W$
- ▲ Min power rating p_{min} (WP) : $3.15W \times 2 = 6.30W$
- ▲ Adjustable angle solar module for ideal sunny condition
- ▲ Physical parameters : Length --- 130.50cm (included aluminum frame)
- ▲ Physical parameters : Width --- 67cm (included aluminum frame)
- ▲ Physical parameters : Depth --- 3.80cm (included aluminum frame)
- ▲ Physical parameters : Weight --- 30kgs (included aluminum frame)
- ▲ Solar power charger : Prostar - 12
- ▲ When batteries were completely dead, it would take approx. 12 hours
- ▲ To fully recharge the batteries under ideal sunny condition
- ▲ Solar powered and AC electric two way use, changeable
- ▲ The other specifications same as power wheelchair model TOP-POWER-A1



TOP-03E Folding Four Wheels Walker

- ▲ Durable aluminum frame
- ▲ Hypalon handgrip and backrest
- ▲ Adjustable handgrip 31"~36"
- ▲ 4" EVA caster
- ▲ Press down brake on rear wheel



TOP-03EA Folding Four Wheels Walker

- ▲ Durable aluminum frame
- ▲ Hypalon handgrip and backrest
- ▲ Adjustable handgrip 31"~36"
- ▲ 5" PU front caster
- ▲ Press down brake on rear wheel



TOP-03EB Folding Four Wheels Walker

- ▲ Durable aluminum frame
- ▲ Hypalon handgrip and backrest
- ▲ Adjustable handgrip 31"~36"
- ▲ 6" PU caster
- ▲ Push to lock brake on rear wheel



TOP-03F Folding Four Wheels Walker

- ▲ Hypalon handgrip and backrest
- ▲ Adjustable handgrip 31"~36"
- ▲ 4" EVA caster
- ▲ Two tip brakes on rear wheel



TOP-03 Folding Delta Walker

- ▲ Light and maneuverable steel walker
- ▲ Hypalon handgrip
- ▲ Adjustable handgrip 30"~36"
- ▲ 8" x 2" pneumatic caster

Options

- With shopping bag (TOP-03+B)
- With three PU 8" x 2" casters



TOP-03+B

TOP-03M Deluxe Folding Walker



- ▲ Large seat with backstrap
- ▲ One touch brake and adjustable height handle
- ▲ 8" x 2" PU caster
- ▲ Adjustable seat height 20"-22"
- ▲ Seat, service tray and shopping basket

TOP-03MA Deluxe Folding Walker



- ▲ Large seat with backstrap
- ▲ One touch brake and adjustable height handle
- ▲ 8" x 2" PU caster
- ▲ Adjustable seat height 20"-22"
- ▲ Seat, service tray and shopping basket

TOP-03DC Child Walker



- ▲ Chrome plated steel adjustable height handgrip
- ▲ Powder coating, folding frame
- ▲ Two 6" solid wheels in front
- ▲ Two rubber tip at rear for braking
- ▲ Adjustable size: (S) 16.50"-21.50"
(M) 18.50"-23.50" (L) 20.50"-25.50"

TOP-03DW Infant Walker



- ▲ Fix sturdy powder coating frame
- ▲ Adjustable height seat
- ▲ Comfort upholstery seat
- ▲ Four 2" casters with brake
- ▲ K/D type for transportation
- ▲ Adjustable : 21"-31"

TOP-0500 Swingdown Arm Commode



- ▲ Heavy duty steel frame
- ▲ Chrome plated or color powder coating
- ▲ Adjustable height : (A)14"-18" (B)17"-21"

TOP-0510 Toilet Assist Commode



- ▲ Feature a raised toilet seat, shower-bath, seat and pail
- ▲ Adjustable height 16"-20"

TOP-0520 Commode With Back



- ▲ Feature a raised toilet seat, shower-bath, seat and pail
- ▲ Full length plastic armrest
- ▲ Adjustable height 16"-20"

TOP-0530 Commode With Caster



- ▲ Heavy duty aluminum frame
- ▲ Full length plastic armrest
- ▲ Four casters, rear two casters with brake
- ▲ Adjustable height 16"-20"

TOP-0540 Folding Commode



- ▲ Heavy duty steel frame with chrome finish
- ▲ Molded plastic armrest
- ▲ Seat width : 18" Seat depth : 20"
- ▲ Adjustable height 16"-20"

TOP-0530A Commode With Back



- ▲ Heavy duty aluminum frame
- ▲ Plastic joint, push button for easy assembly
- ▲ Not tools needed for assembly
- ▲ Adjustable height : (A)16"-20" (B)19"-23"

TOP-0540WA Folding Commode With Caster



- ▲ Heavy duty steel frame with chrome finish
 - ▲ Four casters, rear two casters with brake
 - ▲ Adjustable height 16"-20"
- Options
- Fixed height (TOP-0540W)

TOP-04 Knock Down Shower Chair



- ▲ Aluminum frame
- ▲ Fixed armrest and PVC soft seat
- ▲ Four casters, rear two casters with brake

Adjustable Bath Bench



- ▲ Made of high density polyethylene plastic
 - ▲ Adjustable seat height 13"-17"
- TOP-0550A Bath Bench
TOP-0550B Bath Bench With Back

U-type Adjustable Bath Bench



- ▲ Made of high density polyethylene plastic
 - ▲ Adjustable seat height 13"-17"
- TOP-0550C U-type Bath Bench
TOP-0550D U-type Bath Bench With Back

TOP-0500E Adjustable Shower Chair



- ▲ Made by high density polyethylene plastic
- ▲ Adjustable seat height 13"-17"

Transfer Bath Bench



- ▲ Made by high density
- ▲ Adjustable seat height 13"-17"
- TOP-0500C Transfer Bath Bench
- TOP-0500D Transfer Bath Bench With Back

TOP-0500E Shower Board



- ▲ Anti-slip rubber pad improve the safety
- ▲ Adjustable width : 18"-28"
- Options
- Handgrip

TOP-0660 Toilet Safety Rail



- ▲ Full length plastic armrest
- ▲ Locking bracket attaches with basic toilet seat bolts
- ▲ Adjustable : Height 28"-32" Width 20"-24"

TOP-0650 Raised Toilet Seat



- ▲ Legs made of extruded aluminum
- ▲ Five adjustable support bracket
- ▲ Adjustable height : 3"-6"

TOP-0650A Raised Toilet Seat



- ▲ Made of molded plastic fits on most toilet
- ▲ Attaches easily without tools
- ▲ Height : (A)3" (B)6"

TOP-0650B Raised Toilet Seat With Armrest



- ▲ Made of molded plastic fits on most toilet
- ▲ Fixed armrest and attaches easily without tools
- ▲ Height : (A)3" (B)6"

TOP-0650C Raised Toilet Seat Armrest



- ▲ Made of molded plastic fits on most toilet
- ▲ Flip-up armrest, attache easily without tools
- ▲ Height : (A)4" (B)7"



TOP-0690

Adjustable Bench Stool

- ▲ Removable and height adjustable aluminum legs with suction-tips for easy storage, no tool required while assembling.
 - ▲ Injection plastic round seat is easy clean.
 - ▲ Seat dimension : 32.50cm
- Options
- Round soft seat pad and turning seat with pad provide use comfort



Wall Grab Bar

- ▲ Mount vertical or horizontally for bathing or using commode
- ▲ 1" Steel tube finished with knurled finish for safety
- ▲ Size : 12" ; 18" ; 24"

- TOP-0670** Chrome Plated
- TOP-0670A** Chrome Plated
- TOP-0670B** White Powder Coated



Bath Safety Rail

- ▲ 1" Steel tube finished with knurled finish for safety
- ▲ Easy to installed or removed

- TOP-0680** H 7" x W 18"
- TOP-0680A** H 15" x W 17"
- TOP-0680C** H 8" x W 11"
- TOP-0680B** H 15" Clamp Adjust 3"-7" (White powder coated)



Folding and Adjustable Walker

- ▲ Rugged anodized 1" aluminum shaft with chrome plated
- ▲ Easy to open and close with push button
- ▲ Can be used directly over toilet
- ▲ Adjustable height 32"-36"

TOP-0560A Folding and Adjustable Walker

TOP-0560B Folding and Adjustable Walker

TOP-0560C Folding and Adjustable Walker

Options

- Attach wheel (Refer to OPTIONAL PARTS FOR WALKER)



TOP-0560A



TOP-0560B



TOP-0560C

TOP-0560E Reciprocal Type Folding Walker



- ▲ Same as TOP-0560A but with 7/8" shaft
- ▲ Each side of the frame move alternately giving complete stability, affording maximum freedom of movement

Options

- Attach wheel (Refer to OPTIONAL PARTS FOR WALKER)

TOP-0560H Aluminum Rising Walker



- ▲ With additional lower hand grips for convenient to stand-up from a seated position
- ▲ Adjustable height : (A)31"-38" (B)27"-31"

Options

- Attach wheel (Refer to OPTIONAL PARTS FOR WALKER)

Optional Parts for Walker

T-017 3" wheel ; fit 1" shaft ; H 9"

T-018 3" wheel ; fit 7/8" shaft ; H 11"

T-019 3" wheel ; fit 7/8" shaft ; H 9"

T-020 3" wheel ; fit 7/8" shaft ; H 11"

T-021 3" wheel ; fit 1" shaft ; H 9"

T-022 3" wheel ; fit 1" shaft ; H 11"

T-023 Glide brake ; fit 1" shaft ; H 9"

T-024 Glide brake ; fit 1" shaft ; H 11"

T-025 3" swivel wheel ; fit 1" shaft

T-026 3" swivel wheel ; fit 7/8" shaft

T-027 5" swivel wheel ; fit 1" shaft

T-028 5" swivel wheel ; fit 7/8" shaft



K-Posterior Walker

- ▲ Anodized finish aluminum frame
- ▲ One directional casters with ratchet and allow move mobility
- ▲ Handgrips with flanged end to avoid slip.
- ▲ Adjustable height :
Size : S H 19.50" x W 17.50"
M H 21.50" x W 19.50"
L H 23.50" x W 21.50"
XL H 25.50" x W 23.50"

TOP-05650A With Two Wheels + Two Rubber Stopper

TOP-05650B With Two Wheels + Two Ratch Wheels



Quad Cane

- ▲ Anodized finish aluminum shaft with strong chrome plated steel base
- ▲ Offset handle with hypalon handgrip
- ▲ Easily reverses for right or left hand
- ▲ Adjustable height 29"~38"

TOP-0570A 7" x 11" Base

TOP-0570B 6" x 8" Base

TOP-0570C 5" x 7" Base



Quad Cane

- ▲ Anodized finish aluminum shaft with strong chrome plated steel base
- ▲ Orthopedic handle
- ▲ Easily reverses for right or left hand
- ▲ Adjustable height 29"~38"

TOP-0570D 7" x 11" Base

TOP-0570E 6" x 8" Base

TOP-0570F 5" x 7" Base



Two in One Quad Cane

- ▲ Anodized finish aluminum shaft
- ▲ Heavy duty aluminum base

TOP-0570G With Orthopedic Handle

TOP-0570H With Offset Handle





Pyramid Base Quad Cane

- ▲ Anodized finish aluminum shaft

TOP-0570K Adjustable Height 29"~35"

TOP-0570L Adjustable Height 29"~35"

TOP-0570M Adjustable Height 29"~37"

TOP-0570N Adjustable Height 29"~37"

Steel Pyramid Base Quad Cane

TOP-0570P Adjustable Height 29"~37"



TOP-0610W Aluminum Rising Cane

- ▲ With additional lower hand grips for convenient to stand-up from a seated position
- ▲ Adjustable height : (A)34"~39" (B)31"~36"

Aluminum Rising Cane

- ▲ With additional lower hand grips for convenient to stand-up from a seated position

TOP-0570R Adjustable height

(A)36"~41" (B)33"~38" (C)30"~35"

TOP-0570Q Adjustable height 29"~38"

TOP-0580A Underarm Crutch



- ▲ High strength anodized finish aluminum shaft
- ▲ Rubber armpit pad and handgrip
- ▲ Push button adjusted type
- ▲ Adjustable height : Adult (A)53"~61"
Youth (Y)45"~53"
Child (C)38"~46"

TOP-0580W Wooden Underarm Crutch



- ▲ Made of lauan wood
 - ▲ Rubber armpit pad and handgrip
 - ▲ Adjustable height : Adult (A)49"~55"
Youth (Y)45"~51"
Child (C)35"~42"
- TOP-0580P** Stainless steel bed pan

Forearm Crutch

- ▲ High strength anodized finish aluminum shaft
- ▲ Heavy duty vinyl tip for safety and silence
- ▲ Push button adjusted type
- ▲ Adjustable height : Adult (A)39"~51" Youth (Y)35"~47" Child (C)32"~44"

- TOP-0590A** Metal Cuff and Plastic Silence
TOP-0590B Plastic Cuff and Plastic Silence
TOP-0590C Metal Cuff and Metal Silence



TOP-0590B TOP-0590A TOP-0590C

- ▲ Anatomical handgrip

- TOP-0590D** Plastic Cuff and Metal Silence
TOP-0590E Plastic Cuff and Silence
TOP-0590F Metal Cuff and Silence
TOP-0590G Metal Cuff and Metal Silence



TOP-0590D TOP-0590E TOP-0590F TOP-0590G

Cane and Curutch Attachment

- T-001** 1" ; Rubber tip ; Gray
T-002 7/8" ; Rubber tip ; Gray
T-003 7/8" ; Rubber tip ; Gray
T-004 7/8" ; Rubber tip ; Yellow
T-005 Handgrip ; Yellow
T-006 1" ; Rubber tip ; Gray
T-007 3/4" ; Rubber tip ; Gray
T-008 3/4" ; Rubber tip ; Brown
T-009 Armpit pad ; Yellow
T-010 7/8" ; Rubber tip ; Gray
T-011 1-1/8" ; Rubber tip ; Gray
T-012 5/8" ; Rubber tip ; Gray
T-013 5/8" ; Rubber tip ; Gray
T-014 1/2" ; Rubber tip ; Gray
T-015 1/2" ; Rubber tip ; Brown



TOP-0600 Mini Exerciser



- ▲ Heavy duty rubber leg stopper for safety
- ▲ Adjustable tension control
- ▲ Solid stabilizing frame

Bed Pan



- TOP-0580E** Plastic Bed Pan With Cover
TOP-0580C Plastic Bed Pan
TOP-0580D Plastic Urinal for Man

TOP-0570S Folding Cane Seat



- ▲ Aluminum shaft
- ▲ Easy to carry tripod cane
- ▲ Height : (A) 20" (B) 22"

TOP-043 Split Top Overbed Table



- ▲ Powder coated steel column and base
- ▲ 2" swival caster
- ▲ Adjustable table 24.50"-36"
- ▲ Table size : L 31" x W 16"

TOP-045 Overbed Table With Tilp-Top

- ▲ Powder coated steel column and base
- ▲ 2" swival caster
- ▲ Adjustable table 27"-45"
- ▲ Table size : L 23.50" x W 16"

TOP-044 Automatic Overbed Table

- ▲ Powder coated steel column and base
- ▲ 2" swival caster
- ▲ Adjustable table 27"-41.50"
- ▲ Table size : L 31.50" x W 15"



TOP-045

TOP-044

TOP-004 I.V. Stand



- ▲ Made of stainless steel
- ▲ Adjustable height

TOP-008 Mayo Instrument Tray Trolley



- ▲ Made of stainless steel
- ▲ Adjustable height 800cm-1100cm
- ▲ Tray size 510cm x 380cm x 25cm

TOP-016 Folding Stretcher

- ▲ Made of stainless steel
- ▲ Hard wood handle length 221cm
- ▲ Canvas size : L 182cm x W 60cm



TOP-016A Aluminum Portable Folding Stretcher

- ▲ Lightweight two folds for transportation
- ▲ aluminum poles with rubber handles for improved grip
- ▲ PVC upholstery for ease of cleanliness
- ▲ Aluminum feet prevent the stretcher coming into contact with floor
- ▲ Carry bag protects the stretcher when not in use
- ▲ Size : L 226cm x W 11.50cm x H 56cm



TOP-039 Patient Trolley

- ▲ Stove enamelled steel construction
- ▲ With upholstered cushion, I.V. rod and oxygen holder
- ▲ Four swival casters with brake
- ▲ Size : L 196cm x H 71cm x W 71cm

Options

- With adjustable headrest (TOP-039A)





T-105 Safety Belt for Wheelchair

- (A) Metal Buckle Type
- (B) Plastic Buckle type
- (C) Round Type With Plastic Buckle
- (D) Two Piece with Vecro Adjustable

Special Safety Belt

- ▲ Made by neoprene with nylon knitting for comfort
- ▲ Size : S ; M ; L

T-106 Two End Seat Safety Belt

- ▲ Locking and release by buckle



T-107 Four End Seat Safety Belt

- ▲ Locking and release by buckle



T-108 Star Type Chest Safety Belt

- ▲ Locking and release by four end buckle



T-109 Fully Type Chest Safety Belt

- ▲ Locking and release by buckle



T-110 X Type Chest Safety Belt

- ▲ Locking and release by four end buckle




Secure Deep Back and Comfort Postural Cushion



- (A) Steel hard seatrest panel
- (B) Ventilation nylon cover with pellet gel
- (C) High grade sponge
- (D) Memory ergonomics shape foam

- T-201** Steel Hard Seatrest
- T-202** Seat Memory Foam
- T-203** Steel Hard Backrest
- T-204** Backrest Memory Foam



32



T-220 Adjustable angle and height headrest



Seatrest and Backrest Attachment





THE NAME OF QUALITY

MANY NEW MODELS WILL BE OUT SOON



Top Longmax Products Co., Ltd.

MANUFACTURER & EXPORTER

Office : No. 28, Shih Cheng S. 1st Road, Taichung,
407 Taiwan R.O.C.

P.O. Box : 70-46 Taichung 408 Taiwan R.O.C

Tel : 886-4-2255 -6726 rep 6 lines

Fax : 886-4-2251-9516

e-mail : longmax@ms2.hinet.net

website : <http://www.longmax.com.tw>